



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: SEALING TAPE POLYESTER-STRL CLR Material: polyester acrylate Type: film, gasket, or applicator Quantity per box: 800 Quantity per package: 200 For use with (application): Cell culture, sterile applications Product line: nunctrade; Adhesive: pressure sensitive Certified RNase/DNase free: No Pierceable: No Temperature range (metric): -70 °C to +100 °C Plate Compatibility: All For use with: 250050, 4110-11 Incubation Conditions: General Peelable: Yes Unit Size: Case of 800

Additional Information

No.	FE-0393
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900236366

