



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Type: foil, gasket or applicator Quantity per box: 2 Quantity per package: 2 For use with (application): Handheld applicator for uniform seal Product Line: Nunc™ Adhesive: Pressure Sensitive Certified RNase/DNase Free: No Pierceable: No Plate Compatibility: 96-well plates, 384-well plates Unit Size: Case of 2

Additional Information

No.	FE-0444
Manufacturer (Brand)	Thermo Scientific
VGKL number	290025005

