



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: AI SEALING TAPE FOR 96-WELL PLT Type: foil, gasket or applicator Number of wells: 96 Number per box: 100 Number per package: 100 For use with (application): PCR and storage in 96-well plates Product line: Nunc™ Adhesive: Pressure sensitive Certified RNase/DNase free: Yes Pierceable: Yes Temperature range (metric): -80 °C° to +150 °C° Plate Compatibility: 96-well plates Evaporation Barrier: Moldable, allows visual verification of seal integrity, minimizes evaporation and cross-contamination For use with: 250050, 4110-11 Incubation Conditions: General Peelable: Yes Unit Size: Case of 100

Additional Information

No.	FE-0363
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900232698

