



Thermo Scientific Nunc™ Replication System

Product Images

Description

Nunc Replication System Transfer minute amounts of material from plate to plate or replicate complete DNA libraries with the Thermo Scientific Nunc Replication System. Dimensions (L x W x H): 121.6 x 83.3 x 135.4 cm (4.79 x 3.28 x 5.33 in.) Material: Anodized aluminum with stainless steel needles Type: Nunctrade; Replication System Height (metric): 135.4 mm Length (anglo): 4.79 in. Length (metric): 121.6 mm Height (anglo): 5.33 in. Quantity per carton: 1 Quantity per package: 1 Width (anglo): 3.28 in. Width (metric): 83.3 mm Plate Compatibility: 96-well plate, 384-well plate For use with: 384 wells Format: 384 needles Unit Size: Case of 1

Additional Information

No.	FE-0445
Manufacturer (Brand)	Thermo Scientific
VGKL number	4600250393

