



Thermo Scientific Nunc™ Replication System

Product Images

Description

Nunc Replication System Transfer minute amounts of material from plate to plate or replicate complete DNA libraries with the Thermo Scientific Nunc Replication System. Dimensions (L x W x H): 121.6 x 83.3 x 135.4 cm (4.79 x 3.28 x 5.33 in.) Material: Anodized aluminum with stainless steel needles Type: Nunctrade; Replication System Height (metric): 135.4 mm Length (anglo): 4.79 in. Length (metric): 121.6 mm Height (anglo): 5.33 in. Quantity per carton: 1 Quantity per package: 1 Width (anglo): 3.28 in. Width (metric): 83.3 mm Plate Compatibility: 96-well plate, 384-well plate For use with: 96 wells Format: 96 needles Unit Size: Each

Additional Information

No.	FE-0446
Manufacturer (Brand)	Thermo Scientific
VGKL number	290025052

