



Thermo Scientific Orion™ AQ4500 Turbidimeter

Product Images

Description

Orion™ AQ4500 Turbidimeter Take advantage of the dual LED light source and US EPA test regulatory approved drinking/ wastewater methods offered by the Thermo Scientific™ Orion™ AQ4500 Turbidimeter Light source: LED, multiple Product type: Turbidimeter Auto range selection: Yes Battery life: 2500 h Description: RS232 cable for method uploads and data downloads in AQ4000/AQ4500 Methods with EPA approval: 180.1 Methods with ISO approval: ISO 7027 Certification/Compliance: US EPA approval for Drinking Water Reporting per 40 CFR Part 141, vol. 74, no. 147, Monday, August 3, 2009. Rules and regulations pg 38351-38356.EPA 180.1 and ISO 7030 Unit Size: Each

Additional Information

No.	LD-3254
Manufacturer (Brand)	Thermo Scientific

