



Thermo Scientific Smart-Vue Alarm

Product Images

Description

Smart-Vue Alarm Replaces the use of a USB analog model for issuing Smart-Vue wireless monitoring alarms via voice and/or SMS. Alarm Types: Voice and SMS Display Type: No Display Software Type: SmartVue Server, SmartVue Client Version 2.2 Description: Smart-Vue Alarm, 21-50 sensors Calibration: Not calibrated Product Type: Smart-Vue Alarm Type: Wireless Monitoring Prerequisites: 36-month software license required to send unlimited text and voice call notifications. License based on number of sensors (dual modules count as 2 sensors). Requires SmartVueServer/SmartVue Client version 2.2 or higher. Unit Size: Each

Additional Information

No.	LD-5054
Manufacturer (Brand)	Thermo Scientific

