



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Thermo Scientific Smart-View™ Receiver USB 1.0
or 2.0 Receiver, 25mW+SW 20 sensors - Velcro,
CD**

**€1,530.00
plus VAT &
Shipping**

Product Images



Description

Smart-Vue™ Receivers Make your choice from an assortment of Thermo Scientific™ Smart-Vue™ receivers, including receiver packages with power supply and software for your wireless monitoring solution. Note that you will also need an endpoint module for the monitoring system. Compatible with: Refrigerators, freezers, LN2 storage, CO2 incubators, clean rooms, storage areas, ovens Ordering information: Availability of Smart-Vue products may vary by country. Ask your local Thermo Fisher Scientific representative to provide you with a Smart-Vue wireless monitoring solution tailored to your specific needs. Specifications, terms and conditions, and pricing are subject to change without notice. Not all products are available in all countries. Description: USB 1.0 or 2.0 Receiver, 25mW+SW 20 sensors Includes: Velcro, CD ROM with Smart-Vue client server software 2.0 license up to 20 sensors, software drivers and user manuals. Type: Receiver, 868MHz Certifications/Compliance: Smart-Vue facilitates conformance to 21 CFR Part 11, cGxP and SOP compliance. Operating frequency: 868MHz For use with (devices): Smart-Vuetrade; Wireless Monitoring Solutions Designation: USB 1.0 or 2.0 Receiver, 25mW+SW 20 sensors Includes: Velcro, CD ROM with Smart-Vue client server software 2.0 license up to 20 sensors, software drivers and user manuals. Type: Receiver, 868MHz Certifications/Compliance: Smart-Vue facilitates conformance to 21 CFR Part 11, cGxP and SOP compliance. Operating frequency: 868MHz For use with (devices): Smart-Vuetrade; Wireless Monitoring Solutions

Additional Information

No.	LD-5203
Manufacturer (Brand)	Thermo Scientific

