



Thermo Scientific Sterilin Autoclave Bag

Product Images

Description

Sterilin autoclave bags Sterilin autoclave bags are durable, disposable bags with a blue "biohazard" imprint. These bags are specifically designed for sterilization or disposal of contaminated waste in autoclaves or waste incinerators. Description: Autoclave bags (HDPE) Dimensions (L x W): 810 x 610mm Number per carton: 200 Temperature range (metric): 121deg;C Length (anglo): 32 in. Width (anglo): 24 in. Unit Size: Each

Additional Information

No.	LD-1571
Manufacturer (Brand)	Thermo Scientific

