



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Thermo-Flask™ Liquid Nitrogen Benchtop Container

€519.00
plus VAT &
Shipping

Product Images

Description

Thermo-Flask™ liquid nitrogen benchtop containers Maximize refrigerant retention by using Thermo Scientific™ Thermo-Flask™ liquid nitrogen benchtop containers with borosilicate glass inner vessels. Description: Without lid and handle Capacity (Anglo-American): 0.035 cu. ft. Capacity (metric): 1 l Shipping Weight (Anglo-American): 8 lb. Shipping Weight (metric): 3.6 kg Diameter (metric): 160mm Exterior finish: Stainless steel Inner diameter (metric) Neck: 130mm Type: Tabletop container LNsub2/sub Unit Size: Each

Additional Information

No.	LD-0297
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

