



Thermo Scientific Table Cooling Container

Product Images



Description

Benchtop Cooling Container Protect enzymes, cells, reagents and solutions by maintaining freezing temperatures on the lab bench with the Benchtop Cooling Container. Array: 4 x 5 Construction/Design: Without fill/transparent lid Maintained temp, hours: Hours,Hours,It;-15deg;C, 1 Capacity (metric): 0.2 ml to 0.5 ml Capacity: 20 (0.2 to 0.5 ml tubes) Temperature Range: -20 °C Type: benchtop cooling container Height (metric): 133 mm Height (anglo): 8.25 in. Length (metric): 226 mm Length (anglo): 8.94 in. Width (metric): 165 mm Width (Anglo): 6.5 in. Dimensions (L x W x H): 8.94 x 6.5 x 8.25 in. (226 x 165 x 133 mm) Description: 4x5 chambers,

- 20 degrees Unit Size: Each

Additional Information

No.	FE-0671
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900355501

