



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Trays for Precision Water Baths Test Tube Tray 13-25mm Each Precision

€726.00
plus VAT &
Shipping

Product Images



Description

Trays for Precision Water Baths Hold test tubes and flasks in place with trays for Thermo Scientific™ Precision™ water baths. Test Tube Tray Measures 5 × 10.2in. (12.7 × 25.9cm) Holds 10-13 to 25mm test tubes Flask Trays Measures 12.7 × 25.9cm (5 × 10.2in) 25mL flask tray holds 18 clips 50mL flask tray holds 10 clips Warranty varies by country, please contact us for details Name: Test Tube Tray 13-25mm Type: Tube Tray For use with (devices): Precision Shaking Water Baths Material: Stainless Steel

Additional Information

No.	LD-0503
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

