



Thermo Scientific Reverse Osmosis Membrane

Product Images

Description

Reverse Osmosis Membrane Select the required Thermo Scientific Reverse Osmosis (RO) membrane for your Thermo Scientific Barnstead water treatment system. Description: Reverse Osmosis membrane For use with (units): LabTower EDI 15, LabTower RO 20/60, LabTower TII 20/60, Pacific RO 3/7/12/20, Pacific TII 3/7/12/20 Type: reverse osmosis membrane Unit Size: Each

Additional Information

No.	LD-0344
Manufacturer (Brand)	Thermo Scientific

