



Thermo Scientific Reverse Osmosis Cartridge

Product Images

Description

Reverse Osmosis Cartridge Select the required Thermo Scientific reverse osmosis (RO) cartridge for your Thermo Scientific Barnstead water treatment system. Description: Reverse Osmosis Cartridge, 16 l/h For use with (devices): Smart2Puretrade; Pro Type: Reverse Osmosis Cartridge Unit Size: Each

Additional Information

No.	LD-5502
Manufacturer (Brand)	Thermo Scientific

