



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific VisionMate Highspeed 2D Barcode Reader

€531.00
plus VAT &
Shipping

Product Images



Description

VisionMate Highspeed 2D Barcode Reader The Thermo Scientific VisionMate Highspeed 2D Barcode Reader decodes 2D barcoded supply tubes at high speed – both individually and complete racks. Description: ID- Reader Accessories Quantity: 1/pk, 1/cs. Speed: lt; 1 second Connection type: USB 2.0 port Data transfer method: keyboard wedge Supported barcode types: Code 128, Code 39, Codabar, Interleaved 2 of 5, EAN 13/EAN, EAN 8, EAN-E Type: barcode reader Height (anglo): 5.25 in. Width (anglo): 2.25 in. Depth (anglo): 4 in. Weight (Anglo): 1 lb. For use with (devices): AB-1850 and 312800 VisionMate High Speed Barcode Reader Required Accessories: VisionMate High Speed Barcode Reader and software Unit Size: Each

Additional Information

No.	TA-0467
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

