



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific VisionMate High Speed Barcode Reader

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

VisionMate Highspeed Barcode Reader Quickly read 2D barcodes from tubes in racks of 24, 48, 96 or 384 tubes. The included linear barcode reader helps you track each rack of tubes. Compatible with virtually all database formats. Description: VisionMate Highspeed Barcode Reader Speed: 1 second per rack Data transfer method: CSV, PDF, Excel, ODBC, TCP/IP Decoding capacity: single rack in 24, 48, 96 or 384 tube format ; single tube (in single tube mode) Power supply: 110 to 240 V, 50/60Hz Supported barcode types: Compatible with all 2D barcodes Type: High speed barcode reader Height (metric): 185 mm Height (Anglo): 7.3 inches Width (metric): 130 mm Width (anglo): 5.1 inches Depth (anglo): 7.1 inches Weight (metric): 3.7 kg Length (Anglo): 7.1 inches Length (metric): 180 mm Dimensions (L x W x H): 180 x 130 x 185 mm (7.1 x 5.1 x 7.3 in.) Unit Size: Each

## Additional Information

---

No.	TA-0014
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

