



Thermo Scientific Reusable Cryo Cuff

Product Images

Description

Reusable Cryo Sleeve Use the Thermo Scientific reusable cryo sleeve with tube and canister systems. Description: polyvinyl sleeve; inner diameter x length: 0.6 in (1.5 cm x 10.8 in) (27.4 cm); box of 100 For use with (devices): CryoPlus 1,2,3,4 Type: Cryo Sleeve Material: Polyvinyl Unit Size: Case of 100

Additional Information

No.	LD-0753
Manufacturer (Brand)	Thermo Scientific

