




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific B-Pure™ Water Treatment System

**Price on  
request  
plus VAT &  
Shipping**

### Product Images

---

## Description

---

B-Pure™ Water Treatment System Customize the economical Thermo Scientific™ B-Pure™ laboratory-grade water system to meet the needs of your pretreatment or deionization applications. Certification/Compliance: CE, CSA  
Description: Double holder module with digital purity meter; 240VAC Dimensions (L x W x H): 7 x 15 x 27 in. (17.8 x 38.1 x 68.6 cm) Feed water temperature (Anglo-American): 40 ° to 120 ° F Feedwater temperature (metric): 4 ° to 49 ° C  
Height (Anglo-American): 27 in. Inlet port NPTF (anglo): 0.5 in. Inlet port NPTF (metric): 1.3cm Inlet pressure: 100psig  
Length (Anglo-American): 7 in. Requires: 2 full-size B -Pure cartridges Electrical resistance: 1 to 10MΩ-cm (depending on cartridge used) Type: Double Holder with Digital Purity Meter Width (Anglo-American): 15 in. Flow Rate: 4l/min . Unit Size: Each

## Additional Information

---

No.	LD-3733
Manufacturer (Brand)	Thermo Scientific

