




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific B-Pure™ Water Treatment System

€554.00
plus VAT &
Shipping

Product Images

Description

B-Pure™ Water Treatment System Customize the economical Thermo Scientific™ B-Pure™ laboratory-grade water system to meet the needs of your pretreatment or deionization applications. Certification/Compliance: CE, CSA
Description: Single holder module Dimensions (L x W x H): 17.8 x 17.8 x 61 cm (7 x 7 x 24 in.) Feed water temperature (Anglo-American): 40 °F to 120 °F Feedwater temperature (metric): 4 °C° to 49 °C° Height (Anglo-American): 24 in. inlet port NPTF (Anglo-American): 0.5 in. inlet port NPTF (metric): 1.3 in. Inlet pressure: 100 psig Length (Anglo-American): 7 in. Requires: Full-size B-Pure cartridge Electrical resistance: 1 to 10MΩ-cm (depending on cartridge used)
Type: Single holder water treatment module Width (anglo): 7 in. Flow rate: 4 l/min Unit Size: Each

Additional Information

No.	LD-3731
Manufacturer (Brand)	Thermo Scientific

