



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific TH-641 Swing-out rotor package

**Price on
request
plus VAT &
Shipping**

Product Images

Description

TH-641 Swing-out rotor package Process samples efficiently with the TH-641 swing-out rotor package. Description: TH-641 Swing-out rotor package Type: Swing-out rotor Capacity: 6 x 13.2 ml Max. RZB: 287,660 × g Max. Speed: 41,000 rpm Radius (metric): 153.2 mm Min. radius (metric): 71.9 mm K-factor: 114 Max. Tube length: 89 mm Max. Tube diameter: 14 mm Net weight: 12.3 lb. (5.6 kg) For use with (devices): Thermo Scientific™ Sorvall™ WX ultracentrifuges Angle: 90deg; Average RZB: 211.333 × g Min RZB: 135.005 × g Unit Size: Each

Additional Information

No.	LD-1727
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

