



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific TX-100S clinical swing-out rotor

**€1,880.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

TX-100S clinical swing-out rotor Provide efficient clinical processing of up to 8 x 5/7 ml blood tubes with certified biosafety with these rotors. Description: TX-100S Clinical Swing Out Rotor Type: Swing Out Rotor Capacity: 8 x 15 ml Max. RZB: 3,215 x g Max. Speed: 4,500 rpm Radius (metric): 144 mm K-factor: 14638 Max. Tube length: 106 mm Max. Tube diameter: 17 mm For use with (instruments): Thermo Scientific Sorvall ST 8, Megafuge 8 and SL 8 series centrifuges Angle: 90deg; Biocontainment: Yes Unit Size: Each

## Additional Information

---

No.	LD-2230
Manufacturer (Brand)	Thermo Scientific
VGKL number	107925704
Transport temperature	Room temperature

