



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific TX-750 4 + 750 ml swing-out rotor

€139.00
plus VAT &
Shipping

Product Images

Description

TX-750 4 + 750 ml swing-out rotor Increase laboratory productivity with the Thermo Scientific TX-750 4 x 750 ml swing-out rotor, which offers maximum flexibility at 4700 rpm. Description: Replacement O-rings for lids (round buckets) Type: Replacement O-ring For use with (devices): Thermo Scientific Sorvall Legend XT/XF, Sorvall ST 40, Multifuge X3/X3F, Megafuge 40 and SL 40/40F centrifuges Certification/Compliance: Aerosol tightness certified by CAMR, HPA, Porton Down, United Kingdom Biocontainment: Yes Autoclavable: Not autoclavable Unit Size: Each

Additional Information

No.	LD-2015
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

