



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific Fiberlite 1000mL Bottles Polycarbonate (PC) Pack of 2

**€391.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Fiberlite 1000mL Bottles Thermo Scientific™ 1000mL Fiberlite Bottles combine a Nalgene bottle body with newly designed Fiberlite Cap Closure. These bottles are available for Thermo Scientific™ Sorvall™ LYNX Superspeed and Bios A Large-capacity centrifuges Features - PPCO PC - Increased cycle life compared to previous Thermo Scientific bottles - 100 cycle life (caps + bottles) - Certifications USP Class VI - Non-Cytotoxic - Max speeds match rotors - Ergonomic speed handle design allows much easier bottle handling and removal Designation: Polycarbonate (PC) Packaging: Set of 2 Volume to rim (metric): 1000mL

## Additional Information

---

No.	LD-0005
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

