



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Fiberlite 250mL Bottle Polypropylene (PPCO) Pack of 2 Set of 6

€261.00
plus VAT &
Shipping

Product Images



Description

Fiberlite 250mL Bottle Thermo Scientific™ 250mL Fiberlite Bottles combine a Nalgene bottle body with newly designed Fiberlite Cap Closure. These bottle will be available for Thermo Scientific™ Sorvall™ LYNX Superspeed and General Purpose Centrifuges. Features - PPCO PC - Increased cycle life compared to previous Thermo Scientific bottles - 100 cycle life (caps + bottles) - Certifications USP Class VI - Non-Cytotoxic - Max speeds match rotors - Ergonomic speed handle design allows much easier bottle handling and removal Designation: Polypropylene (PPCO) Packaging: Set of 6 Volume to rim (metric): 250mL

Additional Information

No.	LD-0008
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

