



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® beakers, high form, with graduation and spout, 800 ml

€136.80
plus VAT &
Shipping

Product Images



Description

DURAN® beakers, tall form, with spout, 800 ml Beakers are cylindrical vessels with a straight side wall. Thanks to this shape, beakers are universally applicable laboratory vessels and are usually used as boiling vessels, as collecting vessels and for preparing, mixing and heating solutions. The large opening of the beaker makes it easy to fill. Beakers are available in different designs and shapes: high and low form, with and without spout, with and without scale. Special features and advantages: Precise spout for clean work Easy-to-read scale with line graduation for measuring the contents Burnt-in and therefore permanently durable white printed image Application: heating, diluting solutions, mixing several components, measuring liquid quantities

Additional Information

No.	SD-0225
Manufacturer (Brand)	DURAN
VGKL number	112111653
Transport temperature	Room temperature

