



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Erlenmeyer flask, NS 29/32, 50ml

€200.50
plus VAT &
Shipping

Product Images



Description

DURAN® Erlenmeyer flasks, NS 29/32, 50 ml Due to their large footprint, Erlenmeyer flasks are suitable for preparing and storing solutions, for example. The additional ground joint enables a tight closure of the vessel. Special features and advantages: Tight closure due to standard ground joint Easy-to-read scale with line graduation for measuring the contents Baked-on and therefore permanently durable white printed image The conical shape facilitates mixing Application: Mixing of several components, dilution of solutions, measuring of liquid quantities, heating, Volume sizes 50 ml (NS 24/29), 100 ml (NS 24/29), 200 ml (NS 29/32), 250 ml (NS 45/40), 300 ml (NS 29/32), 500 ml (NS 45/40) and 1000 ml (NS 45/40) are not according to DIN ISO.

Additional Information

No.	SD-1210
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

