



DURAN® Penicillin flask, 4000 ml

Product Images



Description

DURAN® Penicillin flask, 4000 ml Microbiology is an important field of biology and medicine and deals with microorganisms (bacteria, fungi, protozoa or viruses) and their habitats and conditions. For this purpose it is essential to produce pure culture media on/in which the corresponding microorganisms can grow. Irradiation of the media with IR or UV light can accelerate the growth of the cultures. Culture flasks have large bottom surfaces. They should be as flat as possible so that the layer thickness of the medium is uniformly high. Special features and advantages: In this design, the neck is bent upwards, which prevents the liquid culture medium from leaking. The flat, large bottom surface allows the layer thickness of microorganisms to grow evenly high. The large liquid surface ensures better aeration and thus faster growth of microorganisms. Application: preparation of culture media and cultivation, preparation of cultures, liquid culture.

Additional Information

No.	SD-0653
Manufacturer (Brand)	DURAN
VGKL number	192155171
Material	Borosilikatglas DURAN®
Volume	4000 ml

