



Nalgene® Erlenmeyer flask, polycarbonate

Product Images



Description

Nalgene Erlenmeyer Flasks, Polycarbonate Thermo Scientific Nalgene polycarbonate Erlenmeyer flasks are used for the preparation and storage of culture media and for numerous culture procedures. Capacity (anglo): 8 ounces Description: Flask Erlenmeyer Graduated PC 250mL Capacity (metric): 250mL Number Per Case: 12 Number Per Package: 6 Autoclavable: Autoclavable Disposable: No Graduated: Yes Mouth Style: Narrow Base Shape: Plain Cap Size: Stopper No. 6 Unit Size: Case of 12

Additional Information

No.	CD-1419
Manufacturer (Brand)	Nalgene
VGKL number	238741034
Material	Polycarbonat (PC)

