



Nalgene® Erlenmeyer flask, polycarbonate

Product Images



Description

Nalgene Erlenmeyer Flasks, Polycarbonate Thermo Scientific Nalgene polycarbonate Erlenmeyer flasks are used for processing and storing culture media and for numerous culture procedures. Capacity (Anglo): 4 ounces Description: Flask Erlenmeyer Graduated PC 125mL Quantity Per Case: 24 Quantity Per Package: 12 Autoclavable: Autoclavable Disposable: No Graduated: Yes Mouth Style: Narrow Base Shape: Plain Cap Size: Stopper No. 5 Unit Size: Case of 24

Additional Information

No.	CD-1418
Manufacturer (Brand)	Nalgene
VGKL number	238741033
Material	Polycarbonat (PC)

