



## Nalgene® Erlenmeyer flask, polypropylene

### Product Images

---



### Description

---

Nalgene Erlenmeyer Flasks, Polypropylene Thermo Scientific Nalgene polypropylene Erlenmeyer flasks are a shatterproof alternative to glass. Capacity (Anglo): 32 ounces Description: Flask Erlenmeyer PP 1000mL Capacity (metric): 1000mL Material: Polypropylene

Copolymer (PPCO) Quantity Per Case: 6 Quantity Per Package: 2 Autoclavable: Autoclavable Disposable: No Graduated: No Mouth  
Style: Narrow Base Style: Plain Cap Size: Stopper No. 8 Unit Size: Case of 6

## Additional Information

|                      |                   |
|----------------------|-------------------|
| No.                  | CD-1416           |
| Manufacturer (Brand) | Nalgene           |
| Material             | Polypropylen (PP) |

