



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® round bottom flask, narrow neck, with beaded rim, 500 ml

**€135.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

DURAN® round-bottom flasks, narrow neck, 500 ml Bellied vessels are known as flasks in the laboratory sector and, in addition to beakers, are other universally applicable vessels. They are used, for example, for distillation or for heating liquids. Flasks are available in a wide variety of designs: narrow- and wide-necked, with or without ground joint, as round-bottomed or upright flasks. According to the application, the required vessel with its specific components is selected. Special features and advantages: Geometry enables very uniform heating Burnt-in and therefore permanently durable white print Application: use in heating mantles, heating, mixing of several components, sizes 2000, 3000, 5000, 6000, 12000 and 20000 ml not according to DIN ISO, size 12000 ml not according to DIN ISO, according to ASTM E 1403 and size 20000 ml according to ASTM E 1403.

## Additional Information

---

No.	SD-0729
Manufacturer (Brand)	DURAN
VGKL number	112172144
Transport temperature	Room temperature

