



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® round bottom flask, wide neck, with beaded rim, 1000 ml, neck diam. = 50 mm

€25.88 As low
as
€20.96
plus VAT &
Shipping

Product Images



Description

DURAN® round-bottom flasks, wide-neck, 1000 ml Bellied vessels are called flasks in the laboratory sector and, in addition to beakers, are other universally applicable vessels. They are used, for example, for distillation or for heating liquids. Flasks are available in a wide variety of designs: narrow- and wide-necked, with or without ground joint, as round-bottomed or upright flasks. According to the application, the required vessel with its specific components is selected. Special features and advantages: The wide neck enables easy filling and cleaning The geometry enables very uniform heating Burned-in and therefore permanently durable white printed image Application: use in heating mantles, heating, mixing of several components, sizes 50, 1000, 2000, 3000, 5000, 6000, 10000, 12000 and 20000 ml not according to DIN EN ISO.

Additional Information

No.	SD-0750
Manufacturer (Brand)	DURAN
VGKL number	112174154
Transport temperature	Room temperature

