



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Flat bottom flask, narrow neck, with beaded rim, 2000 ml

€47.24 As low
as
€38.25
plus VAT &
Shipping

Product Images



Description

DURAN® Flasks, narrow neck, 2000 ml Bellied vessels are called flasks in the laboratory sector and are, in addition to beakers, other universally applicable vessels. They are used, for example, for distillation or for heating liquids. Flasks are available in a wide variety of designs: narrow- and wide-necked, with or without ground joint, as round-bottomed or upright flasks. According to the application, the required vessel with its specific components is selected. Special features and advantages: Due to the flat bottom, the flask can be set down without a support ring Burnt-in and therefore permanently durable white print Application: heating, sizes 2000, 3000 and 5000 ml not according to DIN ISO.

Additional Information

No.	SD-0720
Manufacturer (Brand)	DURAN
VGKL number	112171164
Transport temperature	Room temperature

