



Bandelin® SONOPULS HD 4050 Homogenizer

Product Images



Description

For volumes from 0.5 - 100 ml Ready-to-use basic equipment for volumes from 0.5 ml to 20 ml consisting of: - Ultrasonic generator GM 4100 - Ultrasonic transducer UW 50 - Titanium sonotrode TS 102, Ø 2 mm SONOPULS homogenizers with amplitude control 10 to 100 % and actual value display, pulsation, AMPLICHRON® circuit (guarantees constant amplitude in independence of changing conditions), integrated timer, remote control. All instruments are CE-marked and are medical devices according to the Directive for In-vitro Diagnostics 98/79/EC. Applications - Digestion of cells, bacteria, viruses, tissue (cell disruption) - Preparation of finest emulsions - Homogenization of substances - Degassing of liquids - Sample preparation for particle size analysis - Acceleration of chemical reactions

Additional Information

No.	BL-0293
Manufacturer (Brand)	Bandelin
VGKL number	100734050

