



FLUICS CONNECT Pay-as-you-go Starter Bundle

Product Images



Description

FLUICS CONNECT Pay-as-you-go

- One workspace license, unlimited users
- Unlimited number of storage objects and locations
- Pairing with up to 2 online label printers RT230

Note: No subscription required, but label packages with pay-as-you-go service for label printing and regular software updates are necessary.

Online Label Printer - Type RT230

- Plug & Play Thermal Transfer Label Printer

- Connectivity: Requires an internet connection to the online print server via LAN cable (RJ45 connector)
- Compatible with all FLUICS online printing solutions, labels, and ribbons
- Warranty: 2 years from the invoice date, except for print heads (6 months) and accessories

Pay-as-you-go Label Package Type 1

- 1,000 labels (1 roll)
- Label specifications: 25.4 x 12.7 mm (1.0 x 0.5"), yellow
- Performance: ultra-permanent, cryogenic, -196°C to +110°C, waterproof, alcohol-resistant
- Requires one roll of R20 ribbon for 4 rolls of labels

Ribbon - Type R20

- 1 roll, 74 m, black
- For cryogenic labels
- Compatible with online label printer type RT230

Pay-as-you-go service with 1,000 labels and regular software updates, guaranteed for at least 2 years from the purchase date.

FLUICS Connect Pro Paket 1 x Drucker

1x PRO Arbeitsbereich-Lizenz - Lizenz auf Lebenszeit - Unbegrenzte Anzahl von Speicherelementen und Speicherorten - Unbegrenzte Anzahl von Benutzern - Mobile und Web-App 1x Online-Etikettendrucker - Plug & Play Thermotransfer-Etikettendrucker - Internetverbindung zum Online-Druckserver über LAN-Kabel (RJ45-Stecker) - 2 Jahre Garantie, ab Rechnungsdatum, ausgenommen Druckköpfe (6 Monate) und Zubehör 1.000 Etiketten Typ 1 - 25,4 x 12,7 mm (1,0 x 0,5"), wasserfest, kryo-beständig von -196°C bis +110°C, ultra-permanent, blanko, gelb - Einschließlich Farbband und Online-Druckservice.

Additional Information

No.	FL-0002
Manufacturer (Brand)	FLUICS
old neoLab article no.	FLCS-FC-PRO-10LP-LT
sterile	No

