



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific P1-CST Multiple Gradient Gel Caster, 10 teeth Standard 10-tooth comb - 1.5mm

€48.40
plus VAT &
Shipping

Product Images



Description

PI-CST Multiple Gradient Gel Caster, 10 teeth Cast one to five gels simultaneously with Fisherbrand™ PI-CST Multiple Gradient Gel Caster, one of the simplest devices to cast gradient gels. With an external gradient maker, gels can be cast with ease and great reproducibility The design and format of this unit will provide years of service in your lab
Designation: Standard 10-tooth comb For use with: PI-CST System Item description: 10 teeth; Tooth thickness, 1.5mm
Number of teeth: 10 Thickness (metric): 1.5mm For use with (devices): Owltrade; PI-CST System Type: Combs

Additional Information

No.	LD-4841
Manufacturer (Brand)	Thermo Scientific

