



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific P2-CST Multiple Gradient Gel Caster Comb, 25 teeth Blank glass

€17.40
plus VAT &
Shipping



Description

P2-CST Multiple Gradient Gel Caster Comb, 25 teeth Cast one to five gels simultaneously with the Fisherbrand™ P2-CST Multiple Gradient Gel Caster, one of the simplest devices to cast gradient gels. With an external gradient maker, gels can be cast with ease and great reproducibility The design and format of this unit will provide years of service in your lab Designation: Blank glass For use with: P10DS-1 Dual Gel Vertical System Item description: Blank glass Thickness of teeth (metric): - Thickness (anglo): 0.18 in. For use with (devices): Fisherbrandtrade; P2-CST System Type: Blank Glass

Additional Information

No.	LD-4868
Manufacturer (Brand)	Thermo Scientific

