



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific P2-CST Multiple Gradient Gel Caster Comb, 25 teeth Standard 15-tooth comb

€48.40
plus VAT &
Shipping

Product Images



Description

P2-CST Multiple Gradient Gel Caster Comb, 25 teeth Cast one to five gels simultaneously with the Fisherbrand™ P2-CST Multiple Gradient Gel Caster, one of the simplest devices to cast gradient gels. With an external gradient maker, gels can be cast with ease and great reproducibility The design and format of this unit will provide years of service in your lab Designation: Standard 15-tooth comb For use with: P2-CST System Item description: Comb, 15 teeth; Tooth thickness, 0.8mm Number of teeth: 15 Thickness of teeth (metric): 0.8mm For use with (devices): Fisherbrandtrade; P2-CST System Type: Combs

Additional Information

No.	LD-4861
Manufacturer (Brand)	Thermo Scientific

