



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific P2-CST Multiple Gradient Gel Caster Comb, 25 teeth, 1.5mm tooth thickness

€48.40
plus VAT &
Shipping

Product Images



Description

P2-CST Multiple Gradient Gel Caster Comb, 25 teeth, 1.5mm tooth thickness Cast one to five gels simultaneously with the Fisherbrand™ P2-CST Multiple Gradient Gel Caster, one of the simplest devices to cast gradient gels. With an external gradient maker, gels can be cast with ease and great reproducibility The design and format of this unit will provide years of service in your lab Designation: Standard 25-tooth comb For use with: P2-CST System Item description: Comb, 25 teeth; Tooth thickness, 1.5mm Number of teeth: 25 Thickness of teeth (metric): 1.5mm For use with (devices): Fisherbrandtrade; P2-CST System Type: Combs

Additional Information

No.	LD-4871
Manufacturer (Brand)	Thermo Scientific

