



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific Owl C2-S Micro Electrophoresis System Horizontal Electrophoresis

**€283.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Owl C2-S Micro Electrophoresis System Screen up to seven samples on one gel in less than five minutes with the Thermo Scientific™ Owl™ C2-S Micro Electrophoresis System. Features: Horizontal microgel system is designed for quick restriction fragment analysis. Microgel may also be used to quickly check a few PCR samples before loading a large number of samples on a larger gel, saving time Glass slides, 2 × 3 in. (5 × 7.5cm), serve as the gel tray To use, pour and form thin gels by surface tension directly on glass slide; complete run lasts five to 10 minutes Using available caster and casting slides, gels can be cast and solidified in less than 10 minutes. Includes: Buffer chamber (buffer volume, 100mL) Lid Power supply leads Casting tray 10 casting slides One seven-well comb, 1.5mm thick Designation: Horizontal Electrophoresis System; Gel size: 5 x 7.5cm Type: Micro Sample range: 7 Buffer volume (metric): 100mL Length Gel (Metric): 7.5cm Width Gel (Anglo): 6.5 Width Gel (Metric): 5cm Height Outside (Anglo): 4.38 in. Length outside (anglo): 6.75 in. Width outside (anglo): 2.63 in. Item description: Horizontal Electrophoresis System; Gel size: 5 x 7.5cm Dimensions (L x W x H) outside: 6.75 x 2.63 x 4.38 in. (17 x 38 x 11cm)

## Additional Information

---

No.	LD-3616
Manufacturer (Brand)	Thermo Scientific

