



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific Owl EasyCast B1 Mini Gel Electrophoresis System

**Price on  
request  
plus VAT &  
Shipping**

### Product Images

## Description

---

Owl EasyCast BI Mini Gel Electrophoresis System Thermo Scientific Owl EasyCast BI mini gel electrophoresis systems are economical enough for any researcher to integrate into their workstation. Gel size: Mini Max. Voltage: 150 V Product Line: Owltrade; Certification/Compliance: CE Marking Description: EasyCast Dimensions (L x W x H): 22 x 15 x 9.5 cm Type: Mini Height (Anglo-American): 3.75 in. Height (metric): 9.5 cm Length (anglo): 8.63 in. Length (metric): 22 cm Length (metric) gel: 11 cm Sample size: 5 to34 Buffer volume (metric): 600 ml Width (Anglo-American): 5.88 in. Width (metric): 15 cm Width (metric) Gel: 9 cm Unit Size: Each

## Additional Information

---

No.	LD-3555
Manufacturer (Brand)	Thermo Scientific

