



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Owl EasyCast B1A Mini Gel Electrophoresis Systems

€435.00
plus VAT &
Shipping

Product Images

Description

Owl EasyCast BIA mini gel electrophoresis systems Thermo Scientific Owl EasyCast BIA mini gel electrophoresis systems provide consistent, trouble-free results. Gel Size: Mini Max. Voltage: 150 V Product Line: Owltrade; Certification/Compliance: CE Marking Description: EasyCast Dimensions (L x W x H): 16 x 10.5 x 9.5 cm Type: Mini Height (Anglo-American): 3.75 in. Height (metric): 9.5 cm Length (anglo): 6.31 in. Length (metric): 16 cm Length (metric) Gel: 8 cm Sample Size: 5 to 24 Buffer Volume (metric): 400 ml Width (Anglo-American): 4.13 in. Width (metric): 10.5 cm Width (metric) Gel: 7 cm Unit Size: Each

Additional Information

No.	LD-3564
Manufacturer (Brand)	Thermo Scientific

