



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Owl EasyCast B3 gel system with built-in recirculation

€680.00
plus VAT &
Shipping

Product Images

Description

Owl EasyCast B3 gel system with built-in recirculation Thermo Scientific Owl EasyCast B3 systems with built-in recirculation are ideal for long sample runs, multiple sample sets or RNA gels. Gel Size: 12 cm x 14 cm Product Line: Owltrade; Certification/Compliance: CE Marking Description: EasyCast; Gel Size: 12 x 14cmL; Buffer volume: 1L; L x W x H: 25 x 18 x 12cm Dimensions (L x W x H): 9.75 x 7 x 4.72 in. (25 x 18 x 12 cm) Type: Mini Height (Anglo-American): 4.72 in. Length (anglo): 9.75 in. Length (metric) Gel: 14cm Sample Size: 8 to 50 Buffer Volume (metric): 1000mL Width (Anglo): 7 in. Width (metric) Gel: 12cm Unit Size: Each

Additional Information

No.	LD-3607
Manufacturer (Brand)	Thermo Scientific

