



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Owl D-2 Wide Gel Electrophoresis Systems

€712.00
plus VAT &
Shipping

Product Images

Description

Owl D-2 Wide Gel Electrophoresis Systems This wide gel range electrophoresis system allows analysis of 10 to 80 samples with only one gel. Max. Voltage: 150 V Product Line: Owltrade; Certification/Compliance: CE Marking Description: complete system with buffer exchange ports Dimensions (L x W x H): 19 x 20 x 10 cm (7.38 x 7.75 x 4 in.) Type: wide system, with buffer exchange ports Height (Anglo-American): 4 in. Height (metric): 10 cm Length (anglo): 7.38 in. Length (metric): 19 cm Length (metric) Gel: 14.4 cm Sample Size: 10 to 80 Buffer Volume (metric): 600 ml Width (Anglo-American): 7.75 in. Width (Anglo-American) Gel: 6.5 Width (metric): 20 cm Width (metric) Gel: 10.2 cm Unit Size: Each

Additional Information

No.	LD-3708
Manufacturer (Brand)	Thermo Scientific

