



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Analytik Jena® Minigel Twin

€1,020.00
plus VAT &
Shipping

Product Images

Description

Minigel-Twin : The Minigel-Twin system is a vertical PAGE electrophoresis system. The robust design provides ease of use and long-term reliability for accurate separation of nucleic acids and proteins. Robust double gel system with thoughtful design. High separation efficiency with short run times. Platinum electrodes and glass plates with fixed glass spacers. Gel size: 8.6 x 7.7 cm. Separation distance: 7.7 cm. Excellent sealing system (profiled silicone seal) ensures user-friendly, leak-proof pouring of polyacrylamide gels. Maximum buffer requirement: 440 ml. Scope of delivery: Basic electrophoresis body including: Buffer chambers, Safety lid, Manual. Note: The system is compatible with: the Eco-Minide polymerized gels, most pre-cast ready-to-use gels (max. 11 x 11 cm). Technical data: Dimensions (width x depth x height): 236 mm x 160 mm x 182 mm; net weight: 2 kg. Software not included.

Additional Information

No.	AJ-0154
Manufacturer (Brand)	Analytik Jena
Transport temperature	Room temperature

