



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Ziegra® CoolNat ice machine, ZBE 150-100 on 40
cm feet capacity 150kg/24h, ice storage 100kg**

**€10,450.00
plus VAT &
Shipping**

Product Images

Description

Ideal ice temperature of $-0.5\text{ }^{\circ}\text{C}$ ZIEGRA-CoolNat ice machine with natural refrigerant,CE compliant, plug-in modular system for connection to Drinking water/VE water (without surcharge) and electricity Standard and micro ice to choose Highest quality and safety through motor protection, water monitoring and automatic shutdown Designed for continuous operation Environmentally friendly technology Hygienic design Comment: Natural refrigerant R290/Propane incl B100 ice storage tank

Additional Information

No.	ZG-0054
Manufacturer (Brand)	Ziegra
Transport temperature	Room temperature

