



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Ziegra® CoolNat ice machine, ZBE 250-200 on
20 cm feet capacity 250kg/24h, ice storage 200
kg**

**€11,270.00
plus VAT &
Shipping**

Product Images

Description

Ideal ice temperature of -0.5°C ZIEGRA-CoolNat ice machine with natural refrigerant, CE compliant, plug-in modular system for connection to Drinking water/VE water (without surcharge) and electricity Standard and micro ice to choose Highest quality and safety through motor protection, water monitoring and automatic shutdown Designed for continuous operation Environmentally friendly technology Hygienic design Comment: Natural refrigerant R290/Propane incl B200 ice storage tank

Additional Information

No.	ZG-0069
Manufacturer (Brand)	Ziegra
Transport temperature	Room temperature

