



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Ziegra® CoolNat ice machine, ZBE 70-35 capacity 70kg/24h, ice storage 35kg

€6,420.00
plus VAT &
Shipping

Product Images

Description

Ideal ice temperature of $-0.5\text{ }^{\circ}\text{C}$ ZIEGRA-CoolNat ice machine with natural refrigerant, CE-compliant, ready to plug in delivery for connection to Drinking water/VE water and electricity Standard and micro ice to choose Highest quality and safety through motor protection, water monitoring and automatic shutdown Designed for continuous operation Environmentally friendly technology Hygienic design Comment: Natural refrigerant R290 / propane

Additional Information

No.	ZG-0045
Manufacturer (Brand)	Ziegra
VGKL number	416607035
Transport temperature	Room temperature

