



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific TSG General Purpose Freezers - Boxed

**€7,330.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

TSG General Purpose Freezers - Boxed The TSG Series of general purpose freezers are built with the performance, security, and quality features you have come to expect from Thermo Scientific. Description: Universal freezer for laboratory applications, TSG Series Temperature range: -12 °C° to -30 °C° Certification/Compliance: CE Capacity (metric): 1500 L Electrical requirements: 230 V/50 Hz Frequency: 50 Hz Plug type: EU Dimensions (L x W x H) external: 837 x 1440 x 2082 mm Dimensions (D x W x H) internal: 690 x 1300 x 1510 mm Shipping weight (metric): 263 kg Refrigerant: R290 Interior: Painted steel Capacity (metric) Rack: 30 kg SmartView Compatibility: Smart-View Pro compatible Warranty: 12 months Temperature (set point): -20 °C Rollers: Yes Connections: 1 Rear cable entry Defrost: Adaptive defrost Control type: Touch screen control with password protection Door design: Fixed Cabinet material: Steel Material: Steel Monitoring options: Alarm No. of Doors: 2 Shelves: 8 adjustable Packaging Type: Packed Unit Size: Each

## Additional Information

---

No.	LD-14992
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

