



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific TSG General Purpose Freezers - Boxed

€5,650.00
plus VAT &
Shipping

Product Images

Description

TSG General Purpose Freezers - Boxed The TSG Series of general purpose freezers are built with the performance, security, and quality features you have come to expect from Thermo Scientific. Description: General purpose freezer for laboratory applications, TSG Series Temperature range: -12 °C° to -30 °C° Certification/Compliance: CE Capacity (metric): 700 L Electrical requirements: 230 V/50 Hz Frequency: 50 Hz Plug type: EU Dimensions (L x W x H) outside: 855 x 750 x 2085 mm Dimensions (D x W x H) inside: 686 x 600 x 1505 mm Shipping weight (metric): 171 kg Refrigerant: R290 Interior: Painted steel Capacity (metric) Shelf: 30 kg Access security: door lock SmartView compatibility: Smart-View Pro compatible Warranty: 12 months Surface treatment exterior: Painted steel Temperature (set point): -20 °C Rollers: Yes Connections: 1 Rear cable entry Defrosting: Adaptive defrost Control type: Digital microprocessor Environmentally friendly features: Energy efficient Door design: Fixed Cabinet material: Steel Material: Steel Monitoring options: Alarm No. of Doors: 1 Shelves: 4 Adjustable Packaging Type: Packed Unit Size: Each

Additional Information

No.	LD-15024
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

